

AVPI1-arginine vasopressin-induced 1 Gene

Reference: PTXBC000877

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PP5395; VIP32; VIT32; arginine vasopressin-induced protein 1; AVP-induced protein 1; vasopressin-induced transcript; arginine vasopressin induced 1		

For research use only.

GENE DESCRIPTION	arginine vasopressin-induced 1
GENE ID	60370
ORIGIN SPECIES	Human
ACCESSIONS	BC000877
VECTOR	pDONR223
NCBI SYMBOL	AVPI1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtacccagcctcgtgtgagcagccccccttgaggcccccattgaggccgggcccgaagcaggcctcagccaacatcttcaggacgcccagctgctgcagatccaaggcctgtttcaacgcagcggggaccagctggccgaggaaacgggcacagatcatctg
ggaaatgtcaggggaccaccgtgtgtgctgagccctcaagagctgcgcaggaagagcccccgaagcagaacccttgggccactgctacaccactgcagccgctcagaatctggagcccccactctgctgccaaccacagagtggccacagagacagctcagtgga
gcagtatctgcactctaggaagaaaagtgccagatcccggaactggaggaagtcaggcccccacaagctacctccaccagatcagacactga