



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

POLR1D-polymerase (RNA) I polypeptide D, 16kDa Gene

Reference: PTXBC000889

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AC19; POLR1C; RPA16; RPA9; RPAC2; RPC16; RPO1-3; TCS2; DNA-directed RNA polymerases I and III subunit RPAC2; DNA-directed RNA polymerase I subunit D; RNA polymerases I and III subunit AC2; polymerase (RNA) I polypeptide D; polymerase (RNA) I polypeptide D		

For research use only.

GENE DESCRIPTION	polymerase (RNA) I polypeptide D, 16kDa
GENE ID	51082
ORIGIN SPECIES	Human
ACCESSIONS	BC000889
VECTOR	pDONR223
NCBI SYMBOL	POLR1D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaaagagagatcaggaagctggagagaaaaataatctggattgaagacctcaatggctgaaggcgagaggaaagacagccctggaatgtctcaggcagctggaaacagatagacactgtgtgacattttattgcacggaagacatacccttagaaattctctacttcatgatcatg
aagaaccggaaagtgaaattttgtgtttacactacgaccattctcagagagcaaaattaatctacgattcaagaatgcaccccttcagctgttgagccattcagagaggccctgaatgagctcatgaattgtctcccaatgtgtcttgacaagttttgaggccagcataaaggactataa
gcatcaaacgcaagtcagaaatccacattctag

More info: Please contact us for any question: contact@proteogenix.fr