

THEM2-thioesterase superfamily member 2 Gene

Reference: PTXBC000894

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	THEM2; HT012; PNAS-27; acyl-coenzyme A thioesterase 13; hypothalamus protein HT012; acyl-CoA thioesterase 13		

For research use only.

GENE DESCRIPTION	thioesterase superfamily member 2
GENE ID	55856
ORIGIN SPECIES	Human
ACCESSIONS	BC000894
VECTOR	pDONR223
NCBI SYMBOL	THEM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgaccagcatgactcagtcctcgtgggaggtgataaaggccatgaccaaggctcgcaattttgagagagttttgggaaagattactctgtctctgctgctcctgggaaagtgattgtgaaatgaaagtagaagaagagcacaatgcaataggcactctccacggcggttgacagcc
acgttagtagataacatacaacaatggctctgctatgcacggaaggaggagcaccggagtcagtgctgatatgaacataacgtacatgtcacctgcacaaattaggagaagatagtgattacagcacatgttctgaagcaaggaaaaaacactgcatttacctctgtggatctgaccaaca
aggccacaggaaaattaatagcacaaggagacacacaaaacacgtgggaaactga