

SSBP1-single-stranded DNA binding protein 1 Gene

Reference: PTXBC000895

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Mt-SSB; SOSS-B1; SSBP; mtSSB; single-stranded DNA-binding protein, mitochondrial; PWP1-interacting protein 17; single-stranded DNA binding protein 1, mitochondrial; single stranded DNA binding protein 1		

For research use only.

GENE DESCRIPTION	single-stranded DNA binding protein 1
GENE ID	6742
ORIGIN SPECIES	Human
ACCESSIONS	BC000895
VECTOR	pDONR223
NCBI SYMBOL	SSBP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtttcgaagacctgtattacaggtacttcgtcagtttgaagacatgagtcggaacaactaccagtttggctcttgaagatccctgaatcgtgtgcacttacttggcgagtggtcaggaccctgtcttgagacaggtggaaggaaaaatccagtcacaatatctctagcaactaatga
gatgtggcgatcaggggatagtgaaagtttaccactgggtgatgtcagtcacaaagacaacatggcacagaatcatcagttatccggcagcctcagagacgtggcatatcaatatgtgaaaaaggggtctcgaatttatttgaagggaataagactatggtgaatacatggataaaataa
tgtgagggcagaagcaacaacatcatagctgataatattatattctgagtgaccagacgaagagaaggagtag