



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

IFITM1-interferon induced transmembrane protein 1 (9-27) Gene

Reference: PTXBC000897

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD225; DSPA2a; IFI17; LEU13; interferon-induced transmembrane protein 1; dispanin subfamily A member 2a; interferon-induced protein 17; interferon-inducible protein 9-27; leu-13 antigen; interferon induced transmembrane protein 1		

For research use only.

GENE DESCRIPTION	interferon induced transmembrane protein 1 (9-27)
GENE ID	8519
ORIGIN SPECIES	Human
ACCESSIONS	BC000897
VECTOR	pDONR223
NCBI SYMBOL	IFITM1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcacaaggagggaacatgagggtggctgctgtgggggcccaccgcacacatctctcaaggtccaccgtgatcaacatccacagcgagacatctcgtgccgaccattctgctgtgctcctgtgtcaaacacctctcttgaactgggtgctgtctgggttcctatagcattgcctactcctgaagcttagggacaggaagatggtgtggcgacgtgaccggggcccgagcctatgcctcaccgcacgaagtgcctgaacatctggggcctgattctgggcatctcatgacgggtgattcatcctgttactctggtattcggctctgtgacagctaccattattgttacagataatcacagaaaaacgggtgtactag

More info: Please contact us for any question: contact@proteogenix.fr