

FLJ22662-hypothetical protein FLJ22662 Gene

Reference: PTXBC000909

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	hypothetical protein FLJ22662
GENE ID	79887
ORIGIN SPECIES	Human
ACCESSIONS	BC000909
VECTOR	pDONR223
NCBI SYMBOL	FLJ22662
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggcagatagtggaagagggtggcagacatctttcaaaatacaactctggcacctatacaatcaatcataggttctggacctgaagaaagtaagctgaaccacagcttgacaaaggcactctgtacattgtggagcaaattctacatatgtagaatattctgaacaaactgatgtt
ctacggaaaggatattggccctcctacaatgttcttccatgaaaaaatctacaactggagtggtatccactgttagttcagaagctgggcttggactactcttatgatttagctccagagccaaatttccggcgtgaccaagggaagtgactgatacggcatccatgaatatatcatgc
gatacaacaattataagaaggatccttacagtagaggtgacccctgtaatacatctgctgctgaggacctgaactcacctaaccgaagtcctggaggttgattgacacaaagggtggcagatatctacatgacatctacacatcctatgccataaagggtgccacagtacaagggtgg
cctccctgttttctgggaccgtttcaacaaactctacatcagggcatggcagaggtctacaacttgattttattaccatgaaccaaatttgaacttgataaaaatga