



Reference: PTXBC000912

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C20orf195; fibronectin type III domain-containing protein 11; fibronectin type III domain containing 11		

For research use only.

GENE DESCRIPTION	chromosome 20 open reading frame 195
GENE ID	79025
ORIGIN SPECIES	Human
ACCESSIONS	BC000912
VECTOR	pDONR223
NCBI SYMBOL	C20orf195
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]