

ZNF611-zinc finger protein 611 Gene

Reference: PTXBC000918

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc finger protein 611		
For research use only.			

GENE DESCRIPTION	zinc finger protein 611
GENE ID	81856
ORIGIN SPECIES	Human
ACCESSIONS	BC000918
VECTOR	pDONR223
NCBI SYMBOL	ZNF611
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgctacaacgattacaatacattggagaatccataatgaagagagatctacaagtgtataaatgtggcaaatgttcagagatcgtcacatattgcaggatcgtggaactcactaggagagaaaccttacaatgtcatgactgtgccagggtcttcagtcaggcttcctatgcaaaa
cataggagaattatacaggagagaaacctcacaagtgtgatgtgtggcaagccttgccttcacgttcacacgtcattagacatcagagaatctatactggacagaaattctacaatgtcatcagtggtggcaagggtcttcagtcagatcactccatgcagaacatcagaaaattcatttt
gagataaccgtcccaatgaaatgagttagcaaacatcaagcattaattgacattggagtcaattcagcattgactga