

SULT1A1-sulfotransferase family, cytosolic, 1A, phenol-preferring, member 1 Gene
Reference: PTXBC000923

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HAST1/HAST2; P-PST; PST; ST1A1; ST1A3; STP; STP1; TSPST1; sulfotransferase 1A1; P-PST 1; aryl sulfotransferase 1; phenol-sulfating phenol sulfotransferase 1; sulfotransferase family, cytosolic, 1A, phenol-preferring, member 1; thermostable phenol sulfotra		

For research use only.

GENE DESCRIPTION	sulfotransferase family, cytosolic, 1A, phenol-preferring, member 1
GENE ID	6817
ORIGIN SPECIES	Human
ACCESSIONS	BC000923
VECTOR	pDONR223
NCBI SYMBOL	SULT1A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagctgatccaggacacctccgccccgacctggagtagctgaaggggtccgctcatcaagtaacttgcagaggcactggggccctgcagagcttcagggccgctgatgacctgtctcatcagcacctaccccaagtcggcactacctgggtaagccagattctggacat
gatctaccaggggtgactggagaaagtgtcacgagctccatctcatcggtgccccttcttgagttcaagccccagggatctccctcagggatggagactctgaaagacacacgggccccagactcctgaagacacacctgcccctggctctgctccccagactctgttgat
cagaaggtcaaggtgtgtatgttcccgcaacgcaaaggtgtggcagtttctactaccactctaccacatggccaaggtgcacccctgagcctgggacacttctctggagaagttcatgttcggagaaggtgtctacggatcctgtaccagcagctgcagaggtgttg
gagctgagccgcaccacacctgttctctaccttctatgaagacatgaaggagaacccgaaagggaattcaaaagatcctggagttgtggggcactccctgcagaggagaccgtggactcgtgtgtcagcacagctgttcaaggagatgaagaagaccctatgaccaactac
accacgtccccaggagttcatgaccacagcatctccccctcatgaggaaggcatggctggggactggaagaccaccttcacgtggcgagaatgagcgttcgatgcggactatgcggagaagatggcaggctgcagcctcagcttccgctctgagctgtga