

AHNAK-AHNAK nucleoprotein Gene

Reference: PTXBC000926

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AHNAK nucleoprotein; AHNAK-related; neuroblast differentiation-associated protein AHNAK; AHNAKRS; PM227; desmoyokin		

For research use only.

GENE DESCRIPTION	AHNAK nucleoprotein
GENE ID	79026
ORIGIN SPECIES	Human
ACCESSIONS	BC000926
VECTOR	pDONR223
NCBI SYMBOL	AHNAK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggagaaggaggagagacaacccgggagctgctgctgccaactggcagggtagtggctccacgggctgaccatgcccagaggagacggcgcttctgtgcaggaggtgacgcagaactccctgcggcccgactgggggtgtcaaggagggggaccagattgtgggt
gccaccatctacttggacaacctgcagtcgggtgaggtgacccagctgctgaacacatggggcaccacacggtggcctgaagctgcaccgcaagggggaccgctctcccgagcctggccagacctggaccgtgaagtcttcagctcctgcagctgaagtggttctgaacacac
cacagccatcagcactggaatgcaagaccagaacaacagaaggaagccagcagccaagccggggcagtttcagctctccaccccaatgcaggactgtag