

TMEM33-transmembrane protein 33 Gene

Reference: PTXBC000948

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	1600019D15Rik; SHINC-3; SHINC3; transmembrane protein 33		
For research use only.			

GENE DESCRIPTION	transmembrane protein 33
GENE ID	55161
ORIGIN SPECIES	Human
ACCESSIONS	BC000948
VECTOR	pDONR223
NCBI SYMBOL	TMEM33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagatagacccgaacggccccaagggcgggcgctgtgcaattcatgatgaccaataaactggacacggcaatgtggctttctcgtttgttcacagttactgtctgtctgtttgttctgcctctcttgggttgcagatgaagcagcaagctttaccaacgtgtttgtggcaaatg
ctcttaccagtgctctgaggtgcatcaagattaccacactccagtttaagcagagcattcctggccaggtttgttagaggacagctgccactacctgtgtattcactcatctttgaaattccatccagttacaatgagatcttcccagctgtttattctcttgcctcatgtgccacatatac
gaaaaaggtccttgacgcaaggggtcaaatagtttacctctgctgagatctgtcttgacaaaataagtgctaatacaaaaaattctgaaattcattgcttgcaatgaaatattcctgatgcctgcgacagttttatgctttttagtggtcaagggaagttgtcccaacctttatatactatagattt
cttacccttcgatattcgtctcgaagaacctatattgtcggaaccttatttaataagactgaggattgtgtgaacacataataatgaaacctgttgcctcctgtttgtgagaagactttgtctccagagcattgcctttataagcagattggcaccacagttccatag