

DYNLT3-dynein, light chain, Tctex-type 3 Gene

Reference: PTXBC000968

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RP3; TCTE1L; CTTEX1L; dynein light chain Tctex-type 3; protein 91/23		

For research use only.

GENE DESCRIPTION	dynein, light chain, Tctex-type 3
GENE ID	6990
ORIGIN SPECIES	Human
ACCESSIONS	BC000968
VECTOR	pDONR223
NCBI SYMBOL	DYNLT3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggagtaccatgccactgcgacgaggttggcttcaatgctgagggaagcccaaatattgtcaagagtggttagatgggttttaggtgtgaagattataatcacaacaacatcaaccagtggactgcaagcatagtggaacaatccttaacacacctggttaagtgggaaagc
ctataaatatattgtgacctgtgcagtggtccagaagagcgcatatggctttcacacagccagctcctgtttttggataccacatctgatgaacctgtaccgtaagatgggagaaccggaccatgaactgtattgtcaacgtttttgccattgctattgttcttaa