

SOX15-SRY (sex determining region Y)-box 15 Gene

Reference: PTXBC000985

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SOX20; SOX26; SOX27; protein SOX-15; SRY (sex determining region Y)-box 15; SRY (sex determining region Y)-box 20; SRY-box 15		

For research use only.

GENE DESCRIPTION	SRY (sex determining region Y)-box 15
GENE ID	6665
ORIGIN SPECIES	Human
ACCESSIONS	BC000985
VECTOR	pDONR223
NCBI SYMBOL	SOX15
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctaccaggtctctacaggaccagcctggagcctggagcctccgggtccacggctgctgctctcatcttcgggacccagcaggcgggagggcgctgggagcccgcgccccgggagcgtgccccgggagaggtgaagcggccgatgaacgcgttcattgtgt
ggagctccgctcagcgcgcagatggcgagcagaacccaagatgcacaactccgagatctccaagcgcctggggcgagtggaagctgctggagcagagcagagcggcccttcgtggaggaggccaagcggctccgccccgacactgcgcactacccgactac
aagtaccggcctcggcgcaaggccaagagctcggcgccggacgttcccgctgcggacagggaagaggcaacctggccagcggcgcccgctctggggcgccgggtacgcgaccaccaaccgagcagagcgtttgggtacagacccccagctactcgcagcctactgcc
tggcagctatggctcttccactgcaaacctggaagccccctaccgtgctcctccctcagagtgacctaggtccagggggaactgctgccacctataccactactgccccctggctcctccactccatacaacctccccctgctgtggtgccccatgccctaaccactctaa