

PRKCH-protein kinase C, eta Gene

Reference: PTXBC001000

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PKC-L; PKCL; PRKCL; nPKC-eta; protein kinase C eta type; protein kinase C eta		

For research use only.

GENE DESCRIPTION	protein kinase C, eta
GENE ID	5583
ORIGIN SPECIES	Human
ACCESSIONS	BC001000
VECTOR	pDONR223
NCBI SYMBOL	PRKCH
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcactagaatagttgataatgaaatgagattttatgaagtataccgtccacatgagcgtctgtctctgtggcctgggatgttaacaggagccaaaaggagggaagtgtgaagaataa