

PDCL3-phosducin-like 3 Gene

Reference: PTXBC001021

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HTPHLP; PHLP2A; PHLP3; VIAF; VIAF1; phosducin-like protein 3; IAP-associated factor VIAF1; VIAF-1; phPL3; viral IAP-associated factor 1; phosducin like 3		

For research use only.

GENE DESCRIPTION	phosducin-like 3
GENE ID	79031
ORIGIN SPECIES	Human
ACCESSIONS	BC001021
VECTOR	pDONR223
NCBI SYMBOL	PDCL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcaggaccccaacgcagacactgaatggaatgacatctacgcaaaagggtatcttaccctcccaaggaaagtctgaagaattggaagaggagcagaagaggagcagcgcacatccacagcagtcagtggtgaaacatatgaagatgactttggaagagctggaggatcatg
aagacgagtttaatgaggaggaatgaacgtgctatgaaatgtacagacggcggagactggctgagtggaagcaactaaactgaagaataaattcggagaagttttggagatctcagggaagattatgttcaagaagtaccaaagctggcggaggctgtgggtcatcttgcacctttac
aaacaagggaattccctctgtgcccgtataaatcagcacctcagtggaacttgcagggaagtttctgatgtcaaatttatcaaaagccatttcaacaacctgcatacccaattatcctgataggaatctgccacgataattgtttaccctggaaggagatatcaaggctcagtttattgtctctgttg
tttggcgcatgacacctgacaagagatgagttggaatggaactgtctgaatctggagcaattatgacagacctggaggaacccctaagaagccgattgaagacgtgttctgtctcagtcggcgctctgtctcatgaagaggacagcgattccgagggtgactga