

RPS10-ribosomal protein S10 Gene

Reference: PTXBC001032

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DBA9; S10; 40S ribosomal protein S10; ribosomal protein S10		
For research use only.			

GENE DESCRIPTION	ribosomal protein S10
GENE ID	6204
ORIGIN SPECIES	Human
ACCESSIONS	BC001032
VECTOR	pDONR223
NCBI SYMBOL	RPS10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgttgatgcctaagaagaaccggattgccattatgaactccttttaaggaggagtcattgtggccaagaaggatgccacatgcctaagcacccggagctggcagacaagaatgtgcccaacctcatgtcatgaaggccatgcagtctcgaagcccgaggctacgtgaaggaca
gtttgccctggagacatttctactgttaccttaccatgagggtatccagtaicctcgtgattaccttcatctgccccggagattgtgcctgccaccctacgccgtagccgtccagagactggcaggcctcgccctaaaggctgtgagggtgagcgacctgcgagactcacaagagggaag
ctgacagagatactacagacggagtgctgtgccacctggtgccgacaagaagccgaggctggggctgggtcagcaaccgaattccagtttagaggcggtttgtgtgtgacgtgtgtcagccacctcagtaa