



Reference: PTXBC001050

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NFAT4; NFATX; nuclear factor of activated T-cells, cytoplasmic 3; NF-ATc3; T cell transcription factor NFAT4; nuclear factor of activated T-cells c3 isoform IE-Xa; nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 3; nuclear factor of activated T-cells, cytoplasmic 3		

GENE DESCRIPTION	nuclear factor of activated T-cells, cytoplasmic, calcineurin-dependent 3
GENE ID	4775
ORIGIN SPECIES	Human
ACCESSIONS	BC001050
VECTOR	pDONR223
NCBI SYMBOL	NFATC3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]