

CLDN7-claudin 7 Gene

Reference: PTXBC001055

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CEPTRL2; CLDN-7; CPETRL2; Hs.84359; claudin-1; claudin-7; clostridium perfringens enterotoxin receptor-like 2; claudin 7		

For research use only.

GENE DESCRIPTION	claudin 7
GENE ID	1366
ORIGIN SPECIES	Human
ACCESSIONS	BC001055
VECTOR	pDONR223
NCBI SYMBOL	CLDN7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaattcgggctgcagttgctgggcttccatggccctgctgggctgggtgggtgctggcctgcacgccatcccgagtgccagatgagctcctatgcgggtgacaacatcatcacggcccaggccatgtacaaggggctgtggatggactgcgtcacgcagagcacgg
ggatgatgagctgcaaaatgtacgactcgggtgctgcccctgctcgcggccttcgaggccactcgaagccctaatgggtgctcctggtgctgggcttcctggccatgtttgtggccacgatgggcatgaagtgcacgcgctgtggggagagacgacaaagtgaagaaggcccgtatagcc
atgggtggaggcataattttcatcgtggcaggtcttgcacacttggtagcttctcctgtatggccatcagattgtcacagactttataaccccttgatccctaccaacattaagtatgagtttggccctgccatcttattggctgggagggctgtccctagtcacctctgggaggtgcactgctct
cctgttctgtcctgggaatgagagcaaggctgggtaccgtgcaccccgctctaccctaagtccaactctccaaggagtatgtga