

HMGB2-high-mobility group box 2 Gene

Reference: PTXBC001063

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	HMG2; high mobility group protein B2; HMG-2; high mobility group protein 2; high-mobility group (nonhistone chromosomal) protein 2; high mobility group box 2		

For research use only.

GENE DESCRIPTION	high-mobility group box 2
GENE ID	3148
ORIGIN SPECIES	Human
ACCESSIONS	BC001063
VECTOR	pDONR223
NCBI SYMBOL	HMGB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtaaaggagacccaacaagccgcgggcaaaatgtcctcgtacgccttcttctgagacgtccgggaagagcacaagaagaacacccggactctcgtcaatttcgcggaatttccaagaagtgttcggagagatggaagaccatgtctgcaaggagaagtgaagt
tgaagatatggcaaaaagtgaacaagctcgtatgacaggagatgaaaaattacgttcctccaaaggtgataagaagggaagaaaaaggacccaatgctctaaaaggccaccatctgccttctctgttttgccttgaaacatgccccaaagatcaaaagtgaacacctgacctat
ccattggggatactgcaagaattgggtgaaatgtggtcagcagtcagccaaagataaacaacatatgaacagaaagcagctaaagctaaaggagaaatgaaaaggatattgctcatatctgccaagggaagtgaaagcagaaagaaggccctggcaggccaacag
gctcaagaagaagaacgaaccagaagatgaggaggaggagggaagaagaagatgaagatgaggagggaaggatgaagatgaagaataa