

MRPL27-mitochondrial ribosomal protein L27 Gene

Reference: PTXBC001066

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	L27mt; 39S ribosomal protein L27, mitochondrial; MRP-L27; mitochondrial ribosomal protein L27		

For research use only.

GENE DESCRIPTION	mitochondrial ribosomal protein L27
GENE ID	51264
ORIGIN SPECIES	Human
ACCESSIONS	BC001066
VECTOR	pDONR223
NCBI SYMBOL	MRPL27
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgtcggtgtgtgtggcgtgaggaccggacagccgttacatccttgaagcccccactccggtacagctcttgcgtcagatacgcacccaagagtcgggtgtagctcctcaaaacctcgggtgaaagtcacaggcagacccaaggcattagaaaatggaagtcactatgttcacgtctgggaacatcattgcaacacagcgccatttccgtggcaccaggtgcccattgtgggtgtgggaagaataatgtctgtatgccctgggaaggggtagtccgctacactaaggaggtctacgtgctcatccagaacacggaggtgtggtatctgaccaggctgccaagggtgctgtctacaagactttgtccacgtgtgtcctgccaagcctgagggcaccttcaactggtagctatgcttga