

C3orf10-chromosome 3 open reading frame 10 Gene

Reference: PTXBC001067

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C3orf10; HSPC300; MDS027; hHBrk1; protein BRICK1; BRICK1, SCAR/WAVE actin-nucleating complex subunit, homolog; haematopoietic stem/progenitor cell protein 300; BRICK1, SCAR/WAVE actin nucleating complex subunit		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 10
GENE ID	55845
ORIGIN SPECIES	Human
ACCESSIONS	BC001067
VECTOR	pDONR223
NCBI SYMBOL	C3orf10
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggacaggagatccggcgagcgggagattcaccaggactgggctaaccgggagtagatgataaccagcagcatcaagaaatcgagactttctcaactcggtcgatatgtctgtcgttcaagactgcaacactaaacgagaaatgacagccctgaacggagaa
tagagtacattgaagctcgggtgacaaaaggtagacactcacctag