

PRNPIP-prion protein interacting protein Gene

Reference: PTXBC001072

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PRNPIP; PINT1; ERI1 exoribonuclease 3; enhanced RNAi three prime mRNA exonuclease homolog 3; prion interactor 1; prion protein-interacting protein; ERI1 exoribonuclease family member 3		

For research use only.

GENE DESCRIPTION	prion protein interacting protein
GENE ID	79033
ORIGIN SPECIES	Human
ACCESSIONS	BC001072
VECTOR	pDONR223
NCBI SYMBOL	PRNPIP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggtggatggtcagccaagcctgcagcaagtgcggagagggtcgaatgatggcgaaggaggcctcttagatccaaagcgaagtcgaatgtttgtcactgtggagactgggactaaaagtcagctcccaggccagtgccagtcacttgggcttgcagtgccgattactcaa
gcagtggaattatctgaaaaggcttacagcttcgcatgggctgctggcccaagaatggactcttagacatgaacaaggcctcagcctgcaacacataggccggccacagcggcattgacgactgcaagaacattgccaacatcatgaagacactgcctatcgaggcttcattc
aagcagacatcgaagccgttctga