

A2LD1-AIG2-like domain 1 Gene

Reference: PTXBC001077

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	A2LD1; gamma-glutamylaminocyclotransferase; AIG2-like domain 1; AIG2-like domain-containing protein 1; gamma-glutamylamine cyclotransferase		

For research use only.

GENE DESCRIPTION	AIG2-like domain 1
GENE ID	87769
ORIGIN SPECIES	Human
ACCESSIONS	BC001077
VECTOR	pDONR223
NCBI SYMBOL	A2LD1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctagcttctgtacggcaccctgaagcgggtcagcccaaccagggctcctcgggagggcggccacggctccgagccttctggcgcgccgcccagcgtggagccctaccgttggtgatcgggggagcacaacatccgtggctgtgcactgcccggtc
ggggcgctctgtggagggcgaggtctacgcggttagacgagcggatgctgccttctggatgacttcgagagttgcccggccctgtaccagcgacggtgctgcgggtacagctgctggaggaccggggcccgccgagaggagccgccagcgccaccggtgcagtgtctt
gtgtacagcaggggccaccttcccgccggagtgggcccgagctcccgacatgacagctacgactccgagggcgccgacggctgcgtacaacccccgggagaacagataa