

C19orf50-chromosome 19 open reading frame 50 Gene

Reference: PTXBC001080

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0459 protein C19orf50; C19orf50; BORCS4; KXDL; MST096; MSTP096; kxDL motif-containing protein 1; KxDL motif containing 1		

For research use only.

GENE DESCRIPTION	chromosome 19 open reading frame 50
GENE ID	79036
ORIGIN SPECIES	Human
ACCESSIONS	BC001080
VECTOR	pDONR223
NCBI SYMBOL	C19orf50
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacctcccgactcggcctcgagggtcttctcgccgcctcctgagcatgtgaacacagatgatgtcaacgccatcatcctggccagaagaacatgctggaccgtttgagaagaccaatgagatgctgctcaactcaacaacctgtccagtgcccgctgcagcagatgagc
gaacgcttctgcaccacacgaggaccctagtagagatgaacgggacctggacagcatctccgcgtatcaggagcgtgaaagggaactggccaggcagcaccagaggccttcagccatccccagaggcatccttctggaggagagatgaagacccatccccag
caccagaccaccattgccacctcagaacagagcagcgggtcatgtgacaccagccccgacacgctctgcctcctcctgagccccggcttcgaggacctgtcccatgtccaggctggctccccagccatcaacggccgagccagacagatgacaggagatgacggcgcaatag