

RASL11B-RAS-like, family 11, member B Gene

Reference: PTXBC001087

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ras-like protein family member 11B; RAS like family 11 member B		
For research use only.			

GENE DESCRIPTION	RAS-like, family 11, member B
GENE ID	65997
ORIGIN SPECIES	Human
ACCESSIONS	BC001087
VECTOR	pDONR223
NCBI SYMBOL	RASL11B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctcattcagaacatgtgaccatcgccgagtagcccgccgggcaacgcccggcctccgactgtgtgtggcgccgcccggcctggtaagatcgccgtggggcgccagcggcggtgggcaagaccgcaactgggtccgggttcctcaccacgattc
gggtgactatgaagaatgcaggtatctctatactagacaagttcagatagaaggtgaacccctggctctcaggttcaagacactccaggtaattcaggtccatgagaacagcctgagctgagtgacagctgaatagggtgattcgtgggcagatgctgtgtgacgttttctccatca
ctgactacaagagctatgaactcatcagccagctccaccagcagtgccagcagctacacctgggcacccggcctggtgtgtgtgtggccaacaaagctgacctgtgtgcacatcaaacaggttgaccctcagcttgacctgagctagccagcatgtaggctgctcattctatgaagt
gtctgtcagtgaaaattataatgatgtctacagcgcttccagctctctgtaaagaggtcagtcacaaacagcagcctagcagtagcaccgagaagcgaagaacctccctcattccaggcccaagtcaccaacatgcaggacctgaagaggaggtttaagcaagccctctgtccaaa
gtgaggactgtcacctccgtctga