

PIH1D1-PIH1 domain containing 1 Gene

Reference: PTXBC001108

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NOP17; PIH1 domain-containing protein 1; nucleolar protein 17 homolog; PIH1 domain containing 1		

For research use only.

GENE DESCRIPTION	PIH1 domain containing 1
GENE ID	55011
ORIGIN SPECIES	Human
ACCESSIONS	BC001108
VECTOR	pDONR223
NCBI SYMBOL	PIH1D1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgaacccgaagctgctgggactggagctaaagcagggcggagggcgatcgggtgctgattcggcgcgattgaggagctgctgctgaggcctcgaaggagctccagcaagcccagacaaccagaccagaatcgacacaaatccagcctcagcctgtgttctgcataaagaccaac
tcctcggaagggaaggttttcatcaacatctgccactccccctctatccctctccccgcgcgacgtgaccgaggaggagctgcttcagatgctagaggaggaccaagctgggtttcgcaccccatgagctctgggagagcctcatgcagaactggatgcaaaaggccagggtgtaccgcc
tacgacgtagctgtcaacagcgacttctaccgggagatgcagaacagcgatttcttgcgggagctcgtgatcacatcgccaggggagggccttgaggacaaatacaacttgcagctgaatccggaatggcgcgatgaagaaccggccatctcgcagcagaacatc
cgctcggagcagcgtcctcgatccaggagctgggggacgtgtacacgccccggggagagctgagtcaggccctgaaaagcctcactgaacctgtggctggaagccccgacctcctctggccgaattgacctcccaactggatggagccctggggctgtcgtggag
atcggggagaaccgctgtgatggggggccccagcagctgtatcatctagacgtttatccccgctgcagatcaactctcatgagagcaaggcagccttccaccggaagagaaagcaattaatggtggccatgccgcttctgctgtgctcttga