

H3F3B-H3 histone, family 3B (H3.3B) Gene

Reference: PTXBC001124

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	H3.3B; histone H3.3; H3 histone, family 3B (H3.3B); H3 histone family member 3B		
For research use only.			

GENE DESCRIPTION	H3 histone, family 3B (H3.3B)
GENE ID	3021
ORIGIN SPECIES	Human
ACCESSIONS	BC001124
VECTOR	pDONR223
NCBI SYMBOL	H3F3B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccgaacaaagcagactgctcgttaagtcacgggtgggaaagcccccgaacagctggccacgaaagccgaggaagagcgtccctctaccggcggggtgaagaagcctcatcgctacagccggggaccgtggcgcttcgagagattcgtcgtttacagaagtcgac
cgagctgctcaccggaagctgccctccagaggtgtgaggagatcgccaggaattcaaaaccgacctgaggtttcagagcgcagccatcggtgcgtgcagagagctagcgaagcgtacctggtgggtctgttcgaagataccaacctgtgtgcatccacgctaagagagtc
ccatcatgcccaagacatccagttggctcgggatacggggagagagagcttaa