



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PPIB-peptidylprolyl isomerase B (cyclophilin B) Gene

Reference: PTXBC001125

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CYP-S1; CYPB; HEL-S-39; OI9; SCYLP; peptidyl-prolyl cis-trans isomerase B; PPIase B; S-cyclophilin; cyclophilin-like protein; epididymis secretory protein Li 39; peptidylprolyl isomerase B (cyclophilin B); rotamase B; peptidylprolyl isomerase B		

For research use only.

GENE DESCRIPTION	peptidylprolyl isomerase B (cyclophilin B)
GENE ID	5479
ORIGIN SPECIES	Human
ACCESSIONS	BC001125
VECTOR	pDONR223
NCBI SYMBOL	PIIB
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr