

DUSP23-dual specificity phosphatase 23 Gene

Reference: PTXBC001140

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DUSP25; LDP-3; MOSP; VH2; dual specificity protein phosphatase 23; VH1-like member Z; VH1-like phosphatase Z; low-molecular-mass dual-specificity phosphatase 3; testicular tissue protein Li 59; dual specificity phosphatase 23		

For research use only.

GENE DESCRIPTION	dual specificity phosphatase 23
GENE ID	54935
ORIGIN SPECIES	Human
ACCESSIONS	BC001140
VECTOR	pDONR223
NCBI SYMBOL	DUSP23
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggtgcagcccccaactctctctgggtgctccggcggtggtgggactggcgctccggcggtccccccactaccagttctgttgacgtggcggtgcggcactggtgtccctgacggagcgcgggccccctcacagcgacagctgccccggcctcacctgcac
cgccctgcgcacccccgactctgcccccgcccccgaccagatcgaccgcttcgtgcagatcgtggacgagcccaacgcacggggagaggctgtggagtgacatgtgctctggcctttggccgactggcaccatgctggcctgttacctggtgaaggagcggggcttggtgca
ggagatgccattgctgaaatccgacgactacgacccggctccatcgagacatgagcaggagaaagcagttctccagttctaccagcgaacgaataa