

SECISBP2-SECIS binding protein 2 Gene

Reference: PTXBC001189

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	selenocysteine insertion sequence-binding protein 2; SECIS binding protein 2		

For research use only.

GENE DESCRIPTION	SECIS binding protein 2
GENE ID	79048
ORIGIN SPECIES	Human
ACCESSIONS	BC001189
VECTOR	pDONR223
NCBI SYMBOL	SECISBP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgctggagggccggggagccgaaagcagggcatcaagttatcagcagatgtcaaacattgtccccagattgccgggctcaatgtggcatggttagagtcctcagaagcatgtgtctccccagctctgcagccacatactatccgtttgtcaggaaaccacagtgacag
aaatgtttactcagtgctggctcccagatctttataaccaaccagttgttaccgaggttttcaaacagtga