



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

HIST3H2A-histone cluster 3, H2a Gene

Reference: PTXBC001193

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	histone H2A type 3; histone 3, H2a; histone cluster 3, H2a; histone cluster 3 H2A		
For research use only.			

GENE DESCRIPTION	histone cluster 3, H2a
GENE ID	92815
ORIGIN SPECIES	Human
ACCESSIONS	BC001193
VECTOR	pDONR223
NCBI SYMBOL	HIST3H2A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctcggctcgttgtaaagcagggttgcaaggcgcgcgccaaggctaaatgtcgcctctgcgcgcggggctgcagttccccctgtggccgcgtgcaccgggttctctcgcgaaggcgaactattcggagcgcgtgggcgccggcccccgtctatctgccgcggtgtctcagtaatt
gactcgcgaatctctggatgttccggcgaacgcgcgcgcgcacacaagaagacgcgcatalctgcgcgcacactgcagctggccatcgcgaacgagagagatctcaacaagctctggggccgcgtgacatctgcgcagggtfctgcctcgtcccaacatcaggccgtactcgtc