

CXorf26-chromosome X open reading frame 26 Gene

Reference: PTXBC001220

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UPF0368 protein Cxorf26; CXorf26; protein PBDC1; 2610029G23Rik; polysaccharide biosynthesis domain-containing protein 1; polysaccharide biosynthesis domain containing 1		

For research use only.

GENE DESCRIPTION	chromosome X open reading frame 26
GENE ID	51260
ORIGIN SPECIES	Human
ACCESSIONS	BC001220
VECTOR	pDONR223
NCBI SYMBOL	CXorf26
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggccaccagtggaaatgatgagccggttccggggagttggtgtctgtggcacatgcgctttctctccagcagagtcgtatggcaacgatcctgacattgagatggcttggccatgagagcaatgcagcatgctgaagtctattacaagctgatttcacagttgaccacagttcc
tgaactcaccaaatagatgacaaatttactctgagttccgaaaaattttgagacccttaggatagatgtgttgaccagaagaactcaagtcagaaatcagccaaagagaagtggaggccattctcttgaagttaattgggattgttgaagactcaactatggtactttgctcgactag
attgttctcagggtacactgaggaaaacacatcttgcggcggatcaattcttgcattgaaattgctcggaaacgggaaggctatacaaaagctgtttatcagttgttcaggacaaagaaggagagaaaaggagtcacaaatggaggagaaaaagagctgacagtggaaga
agagaacaccaagaatggaggagagaaggagctgatgtggagaagaaaaagggaaggaaatcaacagagaagacaaaactgacaaaggaggagaaaaagggaagaagctgacaaagaatcaacaaagtggtaaaagctatgtaa