



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

ZDHHC4-zinc finger, DHHC-type containing 4 Gene

Reference: PTXBC001239

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ZNF374; DHHC-4; zinc finger DHHC domain-containing protein 4; zinc finger protein 374; zinc finger, DHHC domain containing 4; zinc finger DHHC-type containing 4		

For research use only.

GENE DESCRIPTION	zinc finger, DHHC-type containing 4
GENE ID	55146
ORIGIN SPECIES	Human
ACCESSIONS	BC001239
VECTOR	pDONR223
NCBI SYMBOL	ZDHHC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr