

NCBP2-nuclear cap binding protein subunit 2, 20kDa Gene

Reference: PTXBC001255

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CBC2; CBP20; NIP1; PIG55; nuclear cap-binding protein subunit 2; 20 kDa nuclear cap-binding protein; NCBP 20 kDa subunit; NCBP interacting protein 1; cell proliferation-inducing gene 55 protein; nuclear cap binding protein subunit 2, 20kD; nuclear cap bin		

For research use only.

GENE DESCRIPTION	nuclear cap binding protein subunit 2, 20kDa
GENE ID	22916
ORIGIN SPECIES	Human
ACCESSIONS	BC001255
VECTOR	pDONR223
NCBI SYMBOL	NCBP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtcgggtggcctcctgaaggcgtgcgcagcgactcctacgtggagctgagccagctacgggaccagcacttcgggtgacaatgaagaacaagaaaaactgaagaaaagctgtacgttatgttggaaatcttctttacacaactgaagaacaaatctatgaactctcagc
aaaagtgtgtacataaagaaatcattatgggtctggataaatgaagaaacagcatgttgattctgtttgtggaaatattactcagcgcagatgctggaaaacgcatgctgtacataaatggagcgcgtctggtgacccaatcattcgacagactgggacgcagctttaaggaggg
caggcaatacggcgtggcgatctggggccagggttcgggatgagatcggcaggactacgatctgggagaggagctatgaaaactggcagacagaccagtga