



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

MFAP3L-microfibrillar-associated protein 3-like Gene

Reference: PTXBC001279

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NYD-sp9; microfibrillar-associated protein 3-like; microfibrillar-associated protein 3-like; testis development protein NYD-SP9; microfibrillar associated protein 3 like		

For research use only.

GENE DESCRIPTION	microfibrillar-associated protein 3-like
GENE ID	9848
ORIGIN SPECIES	Human
ACCESSIONS	BC001279
VECTOR	pDONR223
NCBI SYMBOL	MFAP3L
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatcgaatgaagagcaatctgactgtgtgcttctacctcttctgcccctttaaactctagatcacctctagaccgcctaagagtggtgactaacagcactttaaaggcactaacgtggctctgggctctgtgcccgttaactgacgaactgaccatacatagctaaagaaaggaacag
tgcctgattaacgtactgtttatgcaacccctgaccacagttcaagtggtatcaatccttgcgaagctgctgaagaagaagaagatgagaaggagagagagaggatgagctttaa