



Reference: PTXBC001288

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD55 molecule (Cromer blood group); CD55 molecule, decay accelerating factor for complement (Cromer blood group); CD55 antigen; CROM; DAF; complement decay-accelerating factor		

For research use only.

GENE DESCRIPTION	CD55 molecule, decay accelerating factor for complement (Cromer blood group)
GENE ID	1604
ORIGIN SPECIES	Human
ACCESSIONS	BC001288
VECTOR	pDONR223
NCBI SYMBOL	CD55
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]