

TMEM109-transmembrane protein 109 Gene

Reference: PTXBC001309

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	transmembrane protein 109; mg23; mitsugumin-23		
For research use only.			

GENE DESCRIPTION	transmembrane protein 109
GENE ID	79073
ORIGIN SPECIES	Human
ACCESSIONS	BC001309
VECTOR	pDONR223
NCBI SYMBOL	TMEM109
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagcctccagcatcattcaccatgggaaagcatgtgtcaaaagccattctgatggctcctagtgcccttctcctccactcagcattggccagtcctgagacatttgaccaccaggccaacagagagagaagcccagttgatgtcttgaccagataggtcgtctg
cgaggacactggatgctggtgattggccagagaccatgacctgggtgacagagctcttctcccaagtgtgtggccatctcatcagccatttctgtgccccttcttctctgtctgggatcgccgcacagctgctgaatgccttgggactagctggtgattacctgcccaggccctgaagc
tcagccctggccaggtccagaccttcctgctgtgggagcagggccctgggtctctactggtgctgtctgtctctcctcggcttggcttggcccttctggtggcgcgacatctgtgggctgaagcttgcattcttccggccgcttctggtccctgatgaggtcgtgctgaccttccac
ccgggccctgctactcctggccttctgctatcctctacgccctgtgagccggtcactggtcctccgagcctctggggcccaactcagggccaaggtgcgagggtggaacgccaggtggaggagctgcctggcgcagaggcgagcgccaagggggcccgaggtgtgaggga
ggagtgga