

ASB8-ankyrin repeat and SOCS box-containing 8 Gene

Reference: PTXBC001321

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PP14212; ankyrin repeat and SOCS box protein 8; ankyrin repeat and SOCS box containing 8		

For research use only.

GENE DESCRIPTION	ankyrin repeat and SOCS box-containing 8
GENE ID	140461
ORIGIN SPECIES	Human
ACCESSIONS	BC001321
VECTOR	pDONR223
NCBI SYMBOL	ASB8
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagttccagatgtgtgtatattatgcagagcattcagagcaataactctctccgagcgttaaccgaacaattgctgccatccgttcctcccatgataatgtagaggacatcagaggggagcagatgtgaactgcactcatggcacactgaagcccttgactgtgcctgtatg
gtgtcagatgctgactgtgtggagttactctggaagagggagcggaggtgaatgccctggatgggtataaccgaacagccctccatgagcagagaaagatgaggtgtgtggaggtcctattggagatgtgtgcaaaccccaatgcttggatggcaacagagataccccacttca
ctgggcagcctttaagaacaatgctgagtggtgtgctgctctagagagcggggcctctgtcaatgccctggattacaacaatgataccgctcagctgggtgctgcatgaaaggaatcttgagagtgctcagcatcctctggattatggcgagaggtcagagtcacacctaatag
gccagacaccatctccgcctgtggctctgtagtcaggggacttggaacagagaaagagactctgtttgagctcctccacagagctgttgacactttgaattgaggaataatggcaccatgccacgagaggtggccagagaccgcagctatgtgaaaaactgactgttctgtg
ctcagctccaggaaactctaaaaacactgctcgtatgctgctgagcctgagcctggagctcagatctccccgatgcagtgaaaggccttcactgccagctcttgaaggaatacctgttacttttagaataag