

FHL3-four and a half LIM domains 3 Gene

Reference: PTXBC001351

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	LIM-only protein FHL3; SLIM2; four and a half LIM domains protein 3; FHL-3; SLIM-2; skeletal muscle LIM-protein 2; four and a half LIM domains 3		

For research use only.

GENE DESCRIPTION	four and a half LIM domains 3
GENE ID	2275
ORIGIN SPECIES	Human
ACCESSIONS	BC001351
VECTOR	pDONR223
NCBI SYMBOL	FHL3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagcagtcattgactgtgcaaatgcaacgagtcctgtatggacgcaagtacatccagacagacagcgccctactgtgtgcccctgctatgacaatactttgccaacacgtgtgctgagtgccagcagcttatcgggcatgactgagggagctgttctatgaagaccgccatttcacagaggctgtctccgctgctgcccgtgccagcgtcactagccgatgaaccttcacctgccaggacagtgagctgctctgcaatgactgctactgcagtgctgttctctcgagtgctccgctgttggggagactgtcatgctgggtcccggaagctggaatggaggccagacatggcatgagcactgcttctgtgagtggtgtgaaacagccactgggctcccgttctttgtgcccagacagggtgctcactactgctgcccctgctatgagaacaagtgtgctctcgtcgtgcccgtgcagcaagacgtgacacaggggtgagtgacataccgtgatcagccgtggcatcgagaatgtctgtgtaccggatgccagacgcccctggcaggcagcagttcacctcccggatgaagatccctactgtgtgcccctgttttggagaactctttgcacctaagtgcagcagctgcaagcggcccatgtaggactcgggtggaggcaagtatgtcctttgaagaccgacactggcaccacaactgttctctcgtgcccgtgctctacccctgtgtgggccagggttcgtaccggatggagaccaagtgtctgcccaggctgtagccagcagggccctaa