

CD151-CD151 molecule (Raph blood group) Gene

Reference: PTXBC001374

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD151 molecule (Raph blood group); hemidesmosomal tetraspanin CD151; CD151 antigen (Raph blood group); CD151 antigen; GP27; MER2; PETA-3; RAPH; SFA1; TSPAN24; membrane glycoprotein SFA-1; platelet surface glycoprotein gp27; platelet-endothelial cell tetra		

For research use only.

GENE DESCRIPTION	CD151 molecule (Raph blood group)
GENE ID	977
ORIGIN SPECIES	Human
ACCESSIONS	BC001374
VECTOR	pDONR223
NCBI SYMBOL	CD151
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtgagttcaacgagaagaagacaacatgtggcaccgtttgcctcaagtacgtgtgtttacctaattgctgtctgtggctggcctggcctgtgtatggcagtgaggcatctggacgtggccctcaagagtgtactacatcagcctgtggcctcaggcaccctacctggccacagcc
tacatcctgtgtgtggcgggcaactgtgtcatgtgtgactgggtctgtggcctgtgtgcccacottcaaggagcgtcggaaactgtgtgacctgtactcctgtctcctcatcatctttctgtggagatcatcgtgtgtatcctgcctacgctactaccagcagctgaacacggagctcaa
ggagaacctgaaggacacatgaccaagcgtaccaccagccgggccaatgagcgtgtgaccagcgtgtggaccagctgcagcaggagtccactgtgtggcagcaacaactcacaggactggcgagacagtgtgtggtatccgtcacaggagccggtggcctgtgtgtccca
gacagctgtcgaagacgggtgtgtgtctttgtggacagcgagaccatgcctccaacatcacaaggtggaggggcgtgtcatcaccaagttggagacattcagcagcagcctgagggctattgggctgtgggagcgcattgctgtgtgcaggtctttgcatgattcacgt
gctgcctgtacaggagttcaagctggagcactactga