



Reference: PTXBC001375

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CN2; CPGL; HEL-S-13; HsT2298; PEPA; cytosolic non-specific dipeptidase; carnosine dipeptidase II; cytosolic nonspecific dipeptidase; epididymis secretory protein Li 13; glutamate carboxypeptidase-like protein 1; peptidase A; CNDP dipeptidase 2 (metallopep		

GENE DESCRIPTION	CNDP dipeptidase 2 (metallopeptidase M20 family)
GENE ID	55748
ORIGIN SPECIES	Human
ACCESSIONS	BC001375
VECTOR	pDONR223
NCBI SYMBOL	CNDP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]