



Reference: PTXBC001382

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ENT1; equilibrative nucleoside transporter 1; equilibrative nitrobenzylmercaptapurine riboside (NBMPR)-sensitive nucleoside transporter; nucleoside transporter, es-type; solute carrier family 29 (equilibrative nucleoside transporter), member 1; solute car		

For research use only.

GENE DESCRIPTION	solute carrier family 29 (nucleoside transporters), member 1
GENE ID	2030
ORIGIN SPECIES	Human
ACCESSIONS	BC001382
VECTOR	pDONR223
NCBI SYMBOL	SLC29A1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]