

PLA2G16-phospholipase A2, group XVI Gene

Reference: PTXBC001387

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AdPLA; H-REV107-1; HRASLS3; HREV107; HREV107-1; HREV107-3; HRSL3; HRAS-like suppressor 3; Ca-independent phospholipase A1/2; H-rev 107 protein homolog; HRAS-like suppressor 1; adipose-specific PLA2; adipose-specific phospholipase A2; group XVI phospholipase		

For research use only.

GENE DESCRIPTION	phospholipase A2, group XVI
GENE ID	11145
ORIGIN SPECIES	Human
ACCESSIONS	BC001387
VECTOR	pDONR223
NCBI SYMBOL	PLA2G16
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcgtgcgccattccagagcctaagcctggagacctgattgagatttttgcctttctacagacactggccatctatgttggcgatgatgatgtgttcctctggccctccaagtgcaggtgcaggagctggtgcagccagtgtcatgtccgcctgactgacaaggccatcgtgaaga
aggaattgctgtatgatgtggccgggagtgacaagtaccaggtcaacaacaacatgatgacaagtactgcgcctgcagcaaaatcatccagcggcgaggagctggtggggcaggaggtgctctacaagctgaccagtgagaactcgagcactttgtgaatgagctgcg
ctatgagctgcggcagtgaccaggtcagagatgtcatcatcgtgcaagcgttgcaagtggtggtgcagccatgagccttattggagtcattgtctcaagaacaagcgacaaaagcaataa