



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PGF-placental growth factor Gene

Reference: PTXBC001422

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

D12S1900; PGFL; PIGF; PLGF; PlGF-2; SHGC-10760; placenta growth factor;

placental growth factor, vascular endothelial growth factor-related protein; placental growth factor

For research use only.

GENE DESCRIPTION	placental growth factor
GENE ID	5228
ORIGIN SPECIES	Human
ACCESSIONS	BC001422
VECTOR	pDONR223
NCBI SYMBOL	PGF
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]