



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

SHISA5-shisa homolog 5 (*Xenopus laevis*) Gene

Reference: PTXBC001463

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

SCOTIN; protein shisa-5; shisa homolog 5; shisa family member 5

For research use only.

GENE DESCRIPTION	shisa homolog 5 (<i>Xenopus laevis</i>)
GENE ID	51246
ORIGIN SPECIES	Human
ACCESSIONS	BC001463
VECTOR	pDONR223
NCBI SYMBOL	SHISA5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggggttcggagcgaacttggccgttggccctgaccatctttgtgctgctgctgactactatcatctgttccacctgctcctgctgctgcttcaagaagctgccgccgaccacgtccgggtgtgacaccaccacatccaccactgtgggtgcatgcccttatctcagctccaaagtgtgc
ccggccagctacctctgagaacagtcatcaggctgctaacacacatgcgctcctcagccaggatgcgcagcagaccctacccaagtgcagtaaccacacccctgacccaggccaccacatgggccaccggcctaccacgagacccctgtgtgaggaggagcagccgcgccctacccgcgcag
ccagctctcttcaaacccgctctatgtagtccctgaagccggccctctga