

TARDBP-TAR DNA binding protein Gene

Reference: PTXBC001487

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ALS10; TDP-43; TAR DNA-binding protein 43; TAR DNA-binding protein-43; TAR DNA binding protein		

For research use only.

GENE DESCRIPTION	TAR DNA binding protein
GENE ID	23435
ORIGIN SPECIES	Human
ACCESSIONS	BC001487
VECTOR	pDONR223
NCBI SYMBOL	TARDBP
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtctgaatatattcgggtaccgaagatgagaacgatgagccattgaataaccatcggaagacgatgggacggtgctgctctccaggttacagccagttccagggcggtgtgggcttcgctacaggaatccagtgctcagtgatgagaggtgtccgctggtagaaggaaattctgcagccccagatgctggctggggaatctgggtatgtgtcaactatccaaaagatacaaaaagaaaaatggatgagacagatgcttcacagcagtgaaagtgaagagcagtcagaaaacatcgatttaatagtgttgggtctccatggaaaacacgaaggacgtgaaagagtatttagtacccttggagaagtcttatggtgcaggtcaagaaagatcttaagactggtcattcaaaaggggtttggctttgtctgtttacggaatatgaacacaaagtgaagtaatgtcacagcgacatgatagatggacgatggtgtgactgcaaaactcctaattctaagcaaaagccaagatgagcctttgagaagcagaaaagtgttggggcgctgtacagaggacatgactgaggatgagctgcccggagtctctctcagtagcgggatgtgatgatgtcttcatcccaagccattcagggcctttgctttgacatttcagatgatcagattgcgcagtcctttgtggagaggacttgatcattaaaggaatcagcgttcataatccaatgcctaa