



Reference: PTXBC001492

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CNTFR-alpha; ciliary neurotrophic factor receptor subunit alpha; CNTF receptor subunit alpha; ciliary neurotrophic factor receptor		

For research use only.

GENE DESCRIPTION	ciliary neurotrophic factor receptor
GENE ID	1271
ORIGIN SPECIES	Human
ACCESSIONS	BC001492
VECTOR	pDONR223
NCBI SYMBOL	CNTFR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]