

CD59-CD59 molecule, complement regulatory protein Gene

Reference: PTXBC001506

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	CD59 molecule, complement regulatory protein
GENE ID	966
ORIGIN SPECIES	Human
ACCESSIONS	BC001506
VECTOR	pDONR223
NCBI SYMBOL	CD59
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggaatccaaggagggtctgtcctgttcgggctgtcgtcgtcgtggtgtcttctgccaatcaggtcatagcctgcagtgctacaactgtcctaaccgaactgctgactgcaaaacagccgtcaattgttcatctgattttgatgcgtgtctcattaccaaagctgggtacaagtgataacaa
gtgttggaagtttgagcattgcaatttcaacgacgtcacaacccgttgagggaataatgagctaactgctgcaagaaggacctgtgtaactttaacgaacagcttgaaaatggtgggacatccttatcagagaaaacagttcttctgctggtgactccatttctggcagcagcctgga
gccttcacccctaa