

F11R-F11 receptor Gene

Reference: PTXBC001533

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD321; JAM; JAM1; JAMA; JCAM; KAT; PAM-1; junctional adhesion molecule A; junctional adhesion molecule 1; platelet F11 receptor; platelet adhesion molecule 1; F11 receptor		

For research use only.

GENE DESCRIPTION	F11 receptor
GENE ID	50848
ORIGIN SPECIES	Human
ACCESSIONS	BC001533
VECTOR	pDONR223
NCBI SYMBOL	F11R
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggacaaaggcgcaagtcgagaggaaactgttgctcttcataatggcgatcctgttgctccctggcattggcgagtggttacagtcacattctgaacctgaagtcagaattcctgagaataatcctgtgaagtgctctgtgctactcggcctttctctcccgtggagtgga
agtttgaccaaggagacaccaccagactcgttgcataataacaagatcacagcttctatgaggaccgggtgacctcttccaactggtatcacctcaagtcctgacacggggaagacactggacatacactgtatggtctctgaggaagcggcaacagctatgggaggtcaag
gtcaagctatcgtctgtgctcctcaagcctacagttaacatccccctctgcccacattgggaaccggcgagtcgtgacatgctcagaacaagatggtccccacctctgaatacactggttcaaatgggatagtgatgcctacgaatcccaaaagcaccctgctcagca
actctctatgctgatcccaacagagagctggtcttggatccccctgacgctctgactggagaatacagctgtgagcagcggatgggtatggacacctgactcaaatcctgtgcatggaagctgtggagcggatgtggggctcatcgtggcagccctctgtaa
ccctgattcctgggaatctggttttggcctcgtggttgcctatagccgagccacttgacagaacaaagaaaggacttcgagtaagaaggtgattacagccagcctagtgcccgaagtgaaggagaatcaaacagacctgcatctctggtgga