

DKK1-dickkopf homolog 1 (*Xenopus laevis*) Gene

Reference: PTXBC001539

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DKK-1; dickkopf-related protein 1; dickkopf 1 homolog; dickkopf-1 like; dickkopf-like protein 1; hDkk-1; dickkopf WNT signaling pathway inhibitor 1		

For research use only.

GENE DESCRIPTION	dickkopf homolog 1 (<i>Xenopus laevis</i>)
GENE ID	22943
ORIGIN SPECIES	Human
ACCESSIONS	BC001539
VECTOR	pDONR223
NCBI SYMBOL	DKK1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatggctctggcgagcgggagctacccgggtcttctgctgagtgaggcgccaccttgaactcgggttcaattcaacgctatcaagaacctgccccaccgctggcgcgctgctggggcaccaggctctgcagt
cagcgccgcgccgggaatcctgtacccggcgagggaataagtagcagaccattgacaactaccagcgtaccgctgcgagagacgaggagtgccgcactgagtagctgcgctagtcccaccgcgagggggacgcaggcgtgcaaatctgtctgcctgcaggagcgccga
aaacgctgcatgctgacgctatgtgctgccccgggaattactgcaaaaatggaatatgtgtgctctgatcaaaatcattccgaggagaaattgaggaaaccatcactgaaagcttggtaatgatcatagcaccttgatgggtattccagaagaaccacctgtctcaaaaatgtatcaca
ccaaaggacaagaaggttctgttctcgggtcatcagactgtgcctcaggattgtgtgtgctagacactctgtccaagatctgtaaacctgtcctgaaagaaggtcaagtgtgtaccaagcataggagaaaaggctctcatggactagaaatattccagcgttgttactgtggagaaggt
ctgtcttccggatagacagaagaatcaccatcaagccagtaattctttaggttcacactgttcagagacactaa